

HOW TO ORDER

META Shock / META Shock X Foot with Spectra Sock

Sandal Toe Foot Shell

Ordering example:
META Shock X Foot, 25 cm, left foot, toe resistance

Møø6 - 25L - 2

Model	Size of Foot	Foot	Toe Resistance
M006 = META Shock X Foot	22-30 cm	L = left R = right	1-9 (Refer to the category chart)

Ordering example:
Sandal Toe Foot Shell, 25 cm, left foot, tan

FSø1 - 25LT

Size of Foot	Side of Foot	Color
22-30 cm	L = left R = right	B = Buff T = Tan M = Medium Brown



The META Lineup

All META feet feature the first-ever unibody foot that combines dynamic response with lightweight, hardware-free design.

META Arc

Polycentric ankle offers unparalleled balance and stability with 20° of side-to-side rotation and 12 mm of shift. For K3-K4 users. PDAC approved for L5981 and L5986.

META Shock

Integrated shock module provides adaptive shock absorption for active K4 users. PDAC approved for L5987.

META Shock X

Integrated shock module with torsional rotation reflexively adapts to movement and reduces torque applied to the residual limb. PDAC approved for L5984 and L5987.

FOR MORE INFORMATION
CALL 800.848.4930 OR 740.869.3377
OR VISIT WILLOWWOOD.COM



Go beyond.
Go META.



Now Available!



It doesn't just push the limits.
It absorbs them.

UNIBODY DESIGN.

Industry-first boltless, unibody carbon fiber design combines light weight with high energy return in an attractive, durable prosthetic.

MAXIMUM DISPLACEMENT. MINIMUM BUILD HEIGHT.

Market-leading levels of deflection from the lowest build height in its class provides optimal agility and patient compatibility.

RESPONSIVE ABSORPTION.

Shock module responds appropriately to every step, delivering vertical shock absorption to reduce impact forces across a broad range of activities.



More degrees of freedom.
Impact protection and rotation.

The META Shock X brings torsion to the high activity META Shock foot, incorporating 40° axial rotation for increased functional mobility and gait fluency with less physical exertion.

TORSIONAL FLEXIBILITY.

Reflexively adapts to variations when walking on uneven terrain and promotes natural rotational movement resulting in less torque applied to the residual limb.



PDAC APPROVED FOR
L-CODES: L5984 and L5987


META SHOCK / META SHOCK X FOOT: SPECIFICATIONS

HEEL HEIGHT: 10 mm (3/8")

APPROXIMATE PRODUCT WEIGHT: 720 g

*Size 26 Category 5 with foot cover and sock

TRIAL PERIOD: 60 days from date of patient fitting



META SHOCK / META SHOCK X FOOT CATEGORY CHART

	Low Impact Level		Medium Impact Level		High Impact Level	
Adjusted User Weight (patient body weight + loads carried)	Sizes 22-24 cm	Sizes 25-30 cm	Sizes 22-24 cm	Sizes 25-30 cm	Sizes 22-24 cm	Sizes 25-30 cm
99-115 lbs (45-51 kg)	1	1	1	1	2	2
116-130 lbs (52-58 kg)	1	1	2	2	3	3
131-150 lbs (59-67 kg)	2	2	3	3	4	4
151-170 lbs (68-76 kg)	3	3	4	4	5	5
171-194 lbs (77-87 kg)	4	4	5	5	-	6
195-221 lbs (88-99 kg)	5	5	-	6	-	7
222-256 lbs (100-115 kg)	-	6	-	7	-	8
257-287 lbs (116-130 kg)	-	7	-	8	-	9
288-324 lbs (131-147 kg)	-	8	-	9	-	-
325-365 lbs (148-165 kg)	-	9	-	-	-	-

