

HOW TO ORDER

META Shock Foot with Spectra Sock

Ordering example:
META Shock Foot, 25 cm, left foot, toe resistance 2, medium impact

Møø5 - 25L - 2M

Sandal Toe Foot Shell

Ordering example:
Sandal Toe Foot Shell, 25 cm, left foot, buff

FSø1 - 25LB

Model	Size of Foot	Foot	Toe Resistance	Impact Level	Size of Foot	Side of Foot	Color
M005 = META Shock Foot	22-28 cm	L = left R = right	1-9 (Refer to the category chart)	L = low M = medium H = high	22-28 cm	L = left R = right	B = Buff T = Tan M = Medium Brown



The META line up

Our META family of products build on the light weight and high-energy return of our first-ever unibody foot design to enhance the lives of users everywhere.

META Arc

Polycentric ankle offers unparalleled range of motion for K3-K4 users. PDAC approved for L5981 and L5986.



META Shock

Integrated shock module provides adaptive shock absorption for active K4 users. PDAC approved for L5987.



FOR MORE INFORMATION
CALL 800.848.4930 OR 740.869.3377
OR VISIT WILLOWWOOD.COM



Go beyond.
Go META.



Now Available!



PDAC APPROVED
LCODE: L5987



It doesn't just push the limits.
It absorbs them.

Introducing META Shock. A high activity foot with a revolutionary unibody design plus an integrated shock module that provides adaptive shock absorption, responding appropriately to any activity. The META Shock's low build height and market-leading deflection-to-build-height ratio provide multiphase shock absorption in a compact, lightweight package.



UNIBODY BREAKTHROUGH.

Industry-first boltless, unibody carbon fiber design combines light weight with high energy return.

MAXIMUM DEFLECTION. MINIMUM BUILD HEIGHT.

Provides class-leading levels of deflection from the lowest build height in its class.

RESPONSIVE ABSORPTION.

Shock module responds appropriately to every step, reducing impact forces across a broad range of activities.


META SHOCK FOOT: SPECIFICATIONS

HEEL HEIGHT: 10 mm (3/8")

APPROXIMATE PRODUCT WEIGHT: 720 g

*Size 26 Category 5 with foot cover and sock

TRIAL PERIOD: 60 days from date of patient fitting



META SHOCK FOOT CATEGORY CHART

Adjusted User Weight (patient body weight + loads carried)	Low Impact Level			Medium Impact Level			High Impact Level		
	Sizes 22-24 cm	Sizes 25-27 cm	Size 28 cm	Sizes 22-24 cm	Sizes 25-27 cm	Size 28 cm	Sizes 22-24 cm	Sizes 25-27 cm	Size 28 cm
99-115 lbs (45-51 kg)	1	1	1	1	1	1	2	2	2
116-130 lbs (52-58 kg)	1	1	1	2	2	2	3	3	3
131-150 lbs (59-67 kg)	2	2	2	3	3	3	4	4	4
151-170 lbs (68-76 kg)	3	3	3	4	4	4	5	5	5
171-194 lbs (77-87 kg)	4	4	4	5	5	5	-	6	6
195-221 lbs (88-99 kg)	5	5	5	-	6	6	-	7	7
222-256 lbs (100-115 kg)	-	6	6	-	7	7	-	8	8
257-287 lbs (116-130 kg)	-	7	7	-	8	8	-	9	9
288-324 lbs (131-147 kg)	-	8	8	-	9	9	-	-	-
325-365 lbs (148-165 kg)	-	9	9	-	-	-	-	-	-

